SPORTSTECH

Akupressur

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Acupressure

In a nutshell!

**Free eBook: Info,
Exercises and Tips for
a healthy body!**

CONTENTS

1. INTRODUCTION: THE ACUPRESSURE MAT - MORE THAN JUST A “GIMMICK”	. 18
2. ORIGIN	. 19
3. WHAT MAKES THE ACUPRESSURE MAT SO SPECIAL?	. 20
4. FOR WHOM IS ACUPRESSURE SUITABLE?	. 21
5. AN EFFECTIVE WAY TO ALLEVIATE MULTIPLE COMPLAINTS	. 22
6. USING YOUR ACUPRESSURE MAT CORRECTLY	. 23
7. STRENGTHENING THE SOUL - THE EFFECT OF THE ACUPRESSURE MAT ON ‘HAPPINESS HORMONES’ (ENDORPHINS) AND NEUROTRANSMITTERS	. 26

1 INTRODUCTION: THE ACUPRESSURE MAT- MORE THAN JUST A “GIMMICK”

Is it really necessary that people with pain, stress, colds or digestive problems have to immediately resort to medication or dietary supplements? It appears there is often no credible alternative and visiting a doctor becomes inevitable. Nevertheless, one thing is clear, the desire and the search for alternative solutions and treatment options is increasing and more and more people are turning to various preventive methods to avoid the use of medication wherever possible.

Physical complaints nowadays are increasingly common and some level of physical discomfort has become part of everyday life. People who work in offices, sit for several hours a day in front of the computer are among those who frequently complain of neck and back complaints, headaches, sleep disturbances and depression.

Along with a desire to be physically fit, for which we often turn to sports and fitness studios, more and more people are turning to and relying on products that they can easily use and that don't require large amounts of their time.

The acupressure mat is an effective alternative to assist with pain relief. The “trend in acupressure mats” not only extends beyond trend setters from Scandinavia, but is now taking the world by storm. This health mat is gaining more and more recognition as a ‘must-have’ lifestyle product and is rapidly becoming almost indispensable.

2 ORIGIN

Acupressure is a well-known and well-tried healing technique originating from the Far East. It works akin to a stimulating massage and is based on a combination of Traditional Chinese Medicine (TCM) and the 'bed-of-nails' of traditional Indian yogis. Although it is closely related to acupuncture, no needles are used at all - rather certain points of the skin of the body are stimulated by pressure alone. Unlike a bed of nails, the western acupressure mat consists of pressure peaks made of high quality ABS plastic, located on needle rosettes. It is actually the same material that Lego bricks are made from. By activating specific pressure points, small stimuli are stimulated which in turn stimulates the body's blood circulation. The acupressure mat is known to stimu-

late the body's circulation thus providing effected surface areas with additional oxygen and essential nutrients which suppress the activity of pain-causing nerve cells. Pain, and indeed other physical ailments, can be relieved by using this method which additionally reduces stress and even supports weight loss. While Acupressure is particularly good for self-help and self-treatment, it does not replace any medical treatment for serious physical complaints.

Important: The acupressure mat is not a 'fix-all' treatment option and in some cases a visit to a doctor remains unavoidable.

3 WHAT MAKES THE ACUPRESSURE MAT SO SPECIAL?

The acupressure mat's scope is as extensive as its effects and an enormous number of physical problem areas can be covered:

1. Self-Treatment

There are, of course, multiple ways to get closer to your ideal health. The acupressure mat, however, is an extremely effective and extremely easy-to-use way to get healthier and feeling better fast. Ultimately, the acupressure mat simply stimulates your body's own internal resources to achieve a really beneficial and positive effect. Exactly how this is achieved is better known by the "scientific" minded, but one thing is clear - acupuncture and acupressure are well-established and effective methods that have been used for millennia to assist the body to achieve optimal balance and relaxation.

But how can the body be so effectively stimulated?

The pressure peaks of the acupressure mats stimulate the nervous system in the body, activate the reflex points and release 'positive feeling' or 'happy' hormones, such as oxytocin. These hormones ensure that stress and pain are alleviated, a level of deep relaxation is achieved and an overall feeling of well-being is established.

2. Relaxation, for example, by meditation.

The relaxing effect of the acupressure mat is evident right from the first time you use it. The relaxing effect is so strong that it is really easy to fall asleep afterwards. For this reason, it is recommended that you use the mat in the evening so that you can end the day feeling optimally well and relaxed. To this end, the acupressure mat can be a significant help in combatting sleep disorders.

3. Sport and Fitness

Even in the domain of competitive sports, the acupressure mat is a perfectly acceptable method to achieve significant performance benefits. An increase in performance is initially accompanied by a certain bodily regeneration phase, which aspires to achieve the optimal physical readiness to be 100% present in competitive situations. Its overuse in sports, however, can cause increasing physical stress. For this reason, using an acupressure mat is the best way to prevent any premature muscle fatigue.

4. Wellness, Beauty, Youthfulness.

Real beauty comes from within and can't be achieved simply by covering over external flaws. In this area too, an immediate connection to relaxation can be made. Those who are most relaxed invariably radiate a sense of inner peace and wellness. If this is enhanced by the release of beneficial substances in the body and the body's blood circulation promoted, it can bring about a radiant appearance as well as a firmer skin.

4 FOR WHOM IS ACUPRESSURE SUITABLE?

There are hardly any restrictions or exceptions when using the acupressure mat. The treatment itself is equally suitable for people of all ages and health. For children, the mat has the same positive effects as for adults and especially at bedtime. As baby skin is more sensitive than adult skin, however, and may be more sensitive to pressure peaks, a cloth of some kind should be placed on top of the mat if necessary. It is important that the acupressure mat is not used as a toy.

For mild complaints, the acupressure mat can be tried as a sole therapy. For more serious illnesses, acupressure may only be used alongwith more conventional medical therapy or treatment and may have a positive supportive effect.

In addition, acupressure can be performed anywhere and without additional aids. It can be enjoyed indoors and out risk-free, without any side effects and is particularly suitable for self-treatment. In case of any severe symptoms, a prompt visit to a doctor is recommended for a thorough check up so that any more serious issues are promptly diagnosed.

5 AN EFFECTIVE WAY TO ALLEVIATE MULTIPLE COMPLAINTS.


Sleep Disorders

Almost every person suffers from insomnia at some time. Factors such as stress, tension and pain can particularly rob us of a decent night's sleep. The area on the body needing treatment lies between the spine and scapula at the heart level and stimulating this area is traditionally used to cure insomnia and emotional tension that impedes proper sleep. The acupressure mat really hits the right spot here and stimulates all the right areas.


All Forms of Depression

Although visiting a doctor is almost inevitable in some circumstances, the acupressure mat has a supportive effect and is beneficial for your mood. Regular and long-term use is recommended so that sufficient well-being can be generated through targeted stimulation and regular release of the benefit-laden endorphins.


Hot flushes

Endorphins can control hot flushes by suppressing the female hormone luteotropin. By stimulating certain stimuli, the released endorphins can calm the nervous system, regulate the heart rhythm and prevent other signs of nervous ailments. Clinically, endorphins are proven to be quite effective body temperature regulators.

Strengthening the Immune System

With no side effects, you can use the acupressure mat to increase the blood flow in the effected body region and can thereby refresh body processes that may be out of balance. Messengers such as 'happiness hormones' are released and flow through the bloodstream. So, you can achieve a maximum effect on your body with a minimum of effort and can increasingly support and improve your health, fitness and vitality.

6 USING YOUR ACUPRESSURE MAT CORRECTLY


Place the mat on a soft surface, on your bed or on a yoga mat on the floor. Each set of Sportstech Acupressure Mats also provides you with an extra cushion that you can use to stimulate the neck and shoulders. Long hair should be pinned up to make the treatment as effective and comfortable as possible.

Lie down on the acupressure mat. Your bare skin should now come into contact with the pressure peaks of the needle rosette. If you find the pressure peaks too strong, then you can put a cloth between the naked skin and the mat to start with or until you becomes accustomed to and comfortable with the sensation of using the mat. Start with a 10-minute daily treatment and increase it by another 10 minutes each day. We recommend a maximum treatment time of 40 minutes per day. Pay attention to how you feel and vary the daily use and application depending on your own judgement; in this way, you can adjust treatment to suit your personal needs and then continue your day refreshed and relaxed.

HOW TO TREAT INDIVIDUAL AREAS OF YOUR BODY AND BENEFIT FROM THE TARGETED USE OF STIMULATION

You do not have to be a professional to use an acupressure mat, but it is always a good idea to use the mat carefully as the pressure peaks can mark the skin and temporarily leave pressure marks.

Lie down with your back on the mat and then breathe deeply and calmly using your abdominal muscles for approx. 10-15 minutes. Afterwards a distinctly pleasant warmth is noticeable and stress slowly and steadily decreases. After approx. 25-40 minutes of relaxation, a sense of well-being follows and this will directly assist you to enjoy a good nights sleep. You have the opportunity to use the acupressure mat for the benefit of your whole body. Since the nerves pass via the back and spinal area, which directly connects to all body parts, they also have a direct influence on the 'feel' of the whole body. You decide what feels comfortable and, therefore, you can decide on which parts of the body you apply the mat to - in order to achieve the best possible treatment. In this way, you target the application of the mat to different body regions, for example, the abdomen, legs, back or buttocks where there are different possibilities and effects to stimulate the metabolism and release or prevent tension.

Back

Lie with your back on the acupressure mat. In order for the lumbar region to come into contact with the pressure peaks, place your heels against the buttocks.

Beneficial Effect: Stimulation of the nerve pathways along the spine

Buttocks

Lie down on the acupressure mat with your buttocks fully on the mat. If you feel too much pressure from the pressure peaks, you can also put a cloth on the mat. If you have already tested it and are comfortable, then you can simply continue the treatment on bare skin.

Beneficial Effect: Relieves tension in the lumbar region and buttock area resulting in more energy and stability. Prevents chronic tension.

SOLES

Put your bare feet on the acupressure mat while sitting in a chair. To increase the pressure, you can also stand on the mat.

Beneficial Effect: The stimulation of the reflex zones on the feet promote blood circulation and has a beneficial influence on different body parts and internal organs.

LEGS

Lay the outstretched legs on the acupressure mat while lying down facing up (supine).

Beneficial Effect: To promote circulation in the legs and prevent varicose veins or similar discomfort.

LEFT AND RIGHT SIDE OF THE FACE

Gently lay your face cheek on the mat. If you feel too much pressure from the pressure peaks, you can place a cloth between the cheek and the acupressure mat.

Beneficial Effect: This treatment causes a release of tension in the jaw region. The resulting increased blood circulation gives your face a fresher complexion and even more radiance and well-being.

ABDOMEN

Lie face down with your stomach on the acupressure mat making sure that the upper edge of the mat is remains under the chest area.

Beneficial Effect: With this process, you release tension in both the diaphragm and respiratory muscles

Treatment with the acupressure mat is particularly suitable during your work breaks. You can increase your efficiency and give your workday fresh impetus. Achieve maximum success and the very best results with regular use of the mat.

Activating numerous acupressure points promotes blood circulation, increases oxygen uptake and better balances the functions of various body parts. 

STRENGTHENING THE SOUL - THE EFFECT OF THE ACUPRESSURE MAT ON 'HAPPINESS HORMONES' (ENDORPHINS) AND NEUROTRANSMITTERS

Treatment using the acupressure mat is proven to stimulate the release of the body's own feel-good hormones, such as endorphins and oxytocins.

These hormones reduce the unwanted effects caused by external stimuli or tension, such as stress and depression.

What exactly happens?


The pressure peaks on the mat located on the needle rosettes stimulate the so-named acupuncture points along the nerves. Since the nerve processes run via the spinal area they affect all parts of the body. This is, therefore, one of the fastest ways to induce a more relaxed state. The impulses pass through the bone marrow into the brain, eventually forming and releasing happiness hormones such as endorphins, oxytocins and other associated, significant substances and hormones in the pituitary gland. This is followed by muscle relaxation, stress reduction and the nervous system working at full speed. After about 10-15 minutes, peptide hormones are released and the intensified level flows through

the body. After another 5 minutes you can feel a significant reduction of pain and stress. Increasing the treatment for up to 25 to 40 minutes makes the pain relief benefits even more noticeable.

Feelings of happiness are from activated nerve messengers that are usually produced only in extreme or exceptional situations. These messengers, however, are also released by using the acupressure mat and have an analgesic, soothing, anxiolytic effect and can even lead to ecstatic states.

If there is a lack of production, lack of release or deficiency of the messenger substances mentioned above, states of depression can easily result.

If you use the acupressure mat daily, the permanent relief from many symptoms will be achieved.

There are numerous scientific studies worldwide that clearly illustrate and attest to the positive health effects of the acupressure mat. Now it is up to us to acknowledge and take advantage of these positive effects and use the mat as an effective aid to help achieve optimal health.

8 THIS MIGHT ALSO INTEREST YOU

AXM400

Cannot get enough of acupressure? Then get the AXM400 with a lining of coconut fiber and an exclusive linen cover. With 8316 pressure peaks, the AXM400 stimulates the blood circulation, supplying any affected area with oxygen and important nutrients. Pain-guiding nerve cells are switched off as if by magic. Use the “vitality mat” for the effective improvement and protection of your well-being and create a blissful, warm and soothing feeling of relaxation.

AXM200

Get the AXM200 and use the acupressure set while lying, sitting or standing and tailor it to your own personal needs. Relieve any tension in your neck, shoulder and back area and stimulate your foot reflex zones. The additional acupressure cushion enables strengthened treatment of the head and neck area.


AXM100

You just can't go wrong with the AXM100 - thanks to the specially designed, innovative hexagonal massage heads and skin-friendly natural materials. Feel full, holistic recovery and find peace in the stress of everyday life. With the acupressure's set of more than 11,000 pressure peaks you can enjoy effective relief from discomfort in the neck, back, buttocks and legs and bring body and mind into balance. With the help of a pressure pad and two cloths, you can increase or decrease the pressure according to your individual needs. The portable bag function enables flexible use on the go - and there are no more limits to your relaxation.

BPX100

This balance pad, made of comfortable and soft durable material, is the ideal 'next level' for your holistic yoga, pilates and fitness training options. Use the Balance Pad BPX100 for a variety of exercises along with coordination training, gymnastics and yoga. This is how to help train your motor skills, your condition, your strength and especially your balance!

**DISCOVER OUR SHOWROOMS IN BERLIN,
DÜSSELDORF OR MUNICH.**

*Test our products at our stores
and buy them on Amazon!*

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Sportstech Brands Holding GmbH