DE

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EN

Acupressure Set Instruction Manual

Scope of delivery

- 1x acupressure mat
- 1x acupressure cushion
- 1x instruction manual
- 1x acupressure towel
- 1x a cupressure back cushion

Application

The acupressure set is very versatile. The most effective way is to use it without clothes from the waist up, so that the back is stimulated particularly strongly.

Beginners should attach the mat to a sofa or armchair back to get used to it and lean carefully with their backrest. In this way you will get a feeling for the acupressure mat while sitting and working on the computer or watching TV.

Advanced users can start using the acupressure mat and cushion while lying on the sofa or bed, for example while reading, meditating or watching TV. This cushions your back and allows you to get used to using the mat without putting too much pressure on your back.

Professionals can place the mat directly on the floor, lie with their back on the mat and distribute the body weight evenly on it. As the hard surface of the floor does not provide a cushioning effect, the back is stimulated more strongly and blood circulation in the area is stimulated. This provides pain relief and relaxation. The recommended duration of use for professionals is 45 minutes and more. This method is particularly suitable for users with chronic pain.

Using the Acupressure Mat and Acupressure Pillow in the evening relaxes the body while reducing stress after a long day and preparing both body and mind for sleep. It is not recommended to sleep on the mat.

Place the acupressure towel on the acupressure mat to reduce the intensity in certain areas of the body. This method is particularly suitable for beginners.

Place the back cushion under the acupressure mat to increase pressure on certain areas of the body. This method is particularly useful when areas such as a hollow back are not reached by the acupressure mat.

Safety instructions

- Keep children and pets away from acupressure mats and pillows.
- The product should not be used with skin irritation, haemophilia or anticoagulants.
- If you are unsure about health problems, always consult a doctor before use.
- Not suitable for children under 5 years of age.
- Use the mat only on a stable surface to avoid injury.

Product care

The product should be washed by hand at 40°. Be careful with the sharp spikes.

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